



## Prior to Use

Soak the electrode for about 10 minutes in tap water to remove dirt and oil stains on the sensor surface. After rinsing the sensor with distilled water, the electrode is ready for use.



## Electrode Maintenance

- Do not touch the platinum black coating on the sensor surface and always keep it clean.



- If there is a build-up of solids inside the sensor, remove carefully, then recalibrate the electrode.
- If you do not use the electrode for long periods, wipe clean with a lint-free tissue and store the electrode in a dry and cool area.
- For the model CON-10, make sure to store the electrode with tap water always.

## Specifications

Model	CON-1	CON-0.1	CON-10
Range	10 $\mu$ S/cm to 20 mS/cm	0 to 100 $\mu$ S/cm	0.1 to 200 mS/cm
Cell Constant	K=1	K=0.1	K=10
Material	Platinum sensor, glass body		
Operating Temperature	0 to 80°C (32 to 176°F)		
Electrode Dimensions	150 (L) $\times$ 12 (Ø) mm (5.9 $\times$ 0.47 in.)		
Cable Length	1 m (3.3 ft.)		
Connector	6-pin mini-DIN		