

P Series Glass pH Electrode

P11

Glass pH electrode, suitable for measuring the general and mildly corrosive solutions



Range	0...14 pH
Operating Temperature	0...80°C/32...176°F
Reference	Ag/AgCl, single junction
Liquid Junction	Annular ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

P11-LiCl

Glass pH electrode, suitable for measuring the non-aqueous solutions



Range	0...14 pH
Operating Temperature	0...80°C/32...176°F
Reference	AgCl, double junction
Liquid Junction	Annular ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

P11-HA

Glass pH electrode, suitable for measuring the high alkalinity solutions



Range	0...14 pH
Operating Temperature	0...80°C/32...176°F
Reference	Ag/AgCl, single junction
Liquid Junction	Annular ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

P13

Glass pH electrode, suitable for measuring the micro-volume samples



Range	0...14 pH
Operating Temperature	0...100°C/32...212°F
Reference	Ag/AgCl, single junction
Liquid Junction	Frit Ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	90(L)×4.3(Ø)mm

P15

Glass pH electrode, suitable for measuring the low conductivity solutions



Range	0...11 pH
Operating Temperature	0...50°C/32...122°F
Reference	Ag/AgCl, double junction
Liquid Junction	Annular ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

P16

Glass pH electrode, suitable for measuring the solutions with Tris buffer



Range	0...14 pH
Operating Temperature	0...80°C/32...176°F
Reference	Ag/AgCl, double junction
Liquid Junction	Frit ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	90(L)×6(Ø)mm

E Series Laboratory pH Electrode

E201-BNC

General purpose pH electrode, suitable for measuring the aqueous solutions



Range	0...14 pH
Operating Temperature	0...60°C/32...140°F
Reference	AgCl, single junction
Liquid Junction	Fiber
Body Type	Epoxy
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

E202-BNC

Flat surface pH electrode, suitable for measuring the semisolids



Range	0...14 pH
Operating Temperature	0...60°C/32...140°F
Reference	AgCl, single junction
Liquid Junction	Teflon
Body Type	Epoxy
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

E203-BNC

3-in-1 pH electrode with built-in temperature sensor, suitable for outdoor applications



Range	0...14 pH
Operating Temperature	0...60°C/32...140°F
Reference	AgCl, single junction
Liquid Junction	Fiber
Body Type	Epoxy
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

E Series Laboratory ORP Electrode

501

Plastic ORP electrode, suitable for measuring the liquids with strong redox potential



Sensor Type	Platinum pin
Operating Temperature	0...60°C/32...140°F
Reference	Ag/AgCl
Liquid Junction	Teflon
Body Type	Epoxy
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

502

Plastic ORP electrode, suitable for measuring the liquids with weak redox potential



Sensor Type	Platinum band
Operating Temperature	0...60°C/32...140°F
Reference	Ag/AgCl
Liquid Junction	Teflon
Body Type	Epoxy
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm

504

Glass ORP electrode, suitable for measuring the mildly corrosive liquids



Sensor Type	Platinum band
Operating Temperature	0...80°C/32...176°F
Reference	Ag/AgCl
Liquid Junction	Annular ceramic
Body Type	Glass
Connector	BNC, 1 m cable
Dimensions	120(L)×12(Ø)mm